

Viral Expression

- ☑ **Overview of Viral Gene Delivery**
- ☑ **AAV (Adeno-Associated Virus) Expression**
- ☑ **Adenoviral Expression**
- ☑ **Lentiviral Expression**
- ☑ **Retroviral Expression**

Recombinant Viral Gene Delivery


Recombinant viral vectors provide a powerful means of delivering a gene into a target cell. There are many viral vectors available, and there are pros and cons to each. Use the following table to select the best viral vector for your research.

Comparison of Viral Vectors for Gene Delivery

	Adeno-Associated Virus (AAV) (p. 41-49)	Adenovirus (p. 50-56)	Lentivirus (HIV-1, FIV, SIV) (p. 57-62)	Retrovirus (MMLV) (p. 63-70)
Gene Expression	Transient or Stable	Transient	Transient or Stable	Stable
Will Infect Dividing Cells	Yes	Yes	Yes	Yes
Will Infect Non-Dividing Cells	Yes	Yes	Yes	No
Integrates into Target Cell Genome	No*	No	Yes	Yes
Immune Response in Target Cells	Very Low	High	Low	Moderate
Relative Viral Titer	XXX	XXXX	XXX	XX
Relative Transduction Efficiency	XXX	XXXX	XXX	XX

*Native AAV can integrate, but recombinant AAV rarely does.

Typical Workflow for Viral Gene Delivery



Cell Biolabs offers kits and reagents for every step in your workflow.

Adeno-Associated Virus Kits & Reagents


Adeno-associated virus (AAV) is less immunogenic than adenovirus or retrovirus. We offer a comprehensive line of AAV kits and reagents to ensure you get the best expression from your AAV expression studies:

- Helper Free Expression Systems
- Helper Free Packaging Systems
- Expression & Control Vectors
- Viral Packaging Cell Line
- Premade AAV Controls
- Purification Kits
- Quantitation / Titer Kit
- Transduction Kits

AAV Helper Free Systems

Production of recombinant AAV requires certain genes from the adenovirus genome, which means that an adenovirus usually needs to be present. The AAV Helper Free System eliminates the need for a helper adenovirus. Most of the required adenoviral genes (E2A, E4 and VA RNA) are provided in a pHelper plasmid, while the required E1 gene is provided by the 293 packaging cells.

- **Safer:** pHelper plasmid eliminates the need for a helper virus
- **Flexible:** Packaging vectors and expression vectors available separately or as one complete system, so you only order what you need
- **Expandable:** All plasmids are provided individually, not in a mixture, so you can amplify in competent cells



Gene Delivery using the AAV Helper Free System.

AAV Helper Free Systems are available for a variety of formats and serotypes:

- **AAV Complete Expression Systems** contain all packaging plasmids plus an expression vector and a GFP control vector: **p. 42-45**
- **AAV Packaging Systems** contain the pHelper plasmid and a serotype-specific Rep-Cap plasmid for use with your own expression construct: **p. 45**
- If you have an AAV packaging system for one serotype and want to try another, choose one of 8 different **AAV Rep-Cap Plasmids** from native serotypes 1 through 6 plus AAV-DJ and AAV-DJ/8: **p. 46**
- If you already have an AAV packaging system and need a cloning vector, choose one of 10 different **AAV Expression Vectors** available individually: **p. 46**
- Want to make a control virus? Choose one of our **AAV Control Vectors**: **p. 47**

AAV Helper Free Complete Expression Systems

AAV Helper Free Complete Expression Systems contain everything you need to produce high-titer recombinant adeno-associated virus:

- pHelper Plasmid
- Rep-Cap Plasmid (serotype specific)
- GFP Control Vector
- Choice of 10 AAV Expression Vectors:
 - Gene Expression (CMV or no promoter)
 - shRNA (U6 promoter)
 - Self complementary (scAAV)

AAV Helper Free Expression Systems are available for the following serotypes:

- Native serotypes 1-6
- AAV-DJ, engineered by DNA family shuffling to form a hybrid capsid from 8 different native serotypes; provides significantly higher infectivity rates *in vitro* (see table below)
- AAV-DJ/8, a mutant of AAV-DJ that exhibits increased uptake in brain and other tissues *in vivo*, similar to serotypes 8 and 9

Cell Line	Cell or Tissue Source	AAV-1	AAV-2	AAV-3	AAV-4	AAV-5	AAV-6	AAV-8	AAV-9	AAV-DJ	AAV-DJ/8
HEK293	Hu Kidney	25	100	2.5	0.1	0.1	5	0.7	0.1	500	0.3
HeLa	Hu Cervix	3	100	2.0	0.1	3.7	1.0	0.2	0.1	667	0.2
HepG2	Hu Liver	3	100	16.7	0.3	1.7	5	0.3	ND	1250	0.5
Hep1A	Ms Liver	20	100	0.2	1.0	0.1	1.0	0.2	0.0	400	0.1
911	Hu Retina	17	100	11.1	0.2	0.1	17	0.1	ND	500	0.0
CHO	Hm Ovary	100	100	14.3	1.4	333	50	10.0	1.0	25000	5.0
COS	Si Kidney	33	100	33	3.3	5.0	14	2.0	0.5	500	0.3
MeWo	Hu Skin	10	100	20	0.3	6.7	10	1.0	0.2	2857	1.0
NIH3T3	Ms Fibroblasts	10	100	2.9	2.9	0.3	10	0.3	ND	500	0.1
A549	Hu Lung	14	100	20	ND	0.5	10	0.5	0.1	1000	0.1
HT1180	Hu Fibroblasts	20	100	10.0	0.1	0.3	33	0.5	0.1	333	0.2

Relative Infectivity Rates of AAV Serotypes. Normalized to AAV-2 = 100. ND = Not determined.

Recent Product Citations

1. Altemeier, W.A. et al. (2012). Transmembrane and extracellular domains of syndecan-1 have distinct functions in regulating lung epithelial migration and adhesion. *J. Biol. Chem.* **287**:34927-34935. (VPK-410-DJ)
2. Zhao, H. et al. (2014). SCAMP5 plays a critical role in synaptic vesicle endocytosis during high neuronal activity. *J. Neurosci.* **34**:10085-10095. (VPK-410-SER1)
3. Moshiri, F. et al. (2014). Inhibiting the oncogenic mir-221 by microRNA sponge: toward microRNA-based therapeutics for hepatocellular carcinoma. *Gastroenterol. Hepatol. Bed Bench* **7**:43-54. (VPK-418-DJ)

AAV-DJ Helper Free Complete Expression Systems

Product Name	Size	Catalog Number
AAV-DJ Helper Free Expression System	1 kit	VPK-410-DJ
AAV-DJ Helper Free Bicistronic Expression System (Puro)	1 kit	VPK-415-DJ
AAV-DJ Helper Free Bicistronic Expression System (Neo)	1 kit	VPK-416-DJ
AAV-DJ Helper Free Bicistronic Expression System (Hygro)	1 kit	VPK-417-DJ
AAV-DJ Helper Free Bicistronic Expression System (GFP)	1 kit	VPK-418-DJ
AAV-DJ Helper Free Bicistronic Expression System (Blasticidin)	1 kit	VPK-419-DJ
AAV-DJ Helper Free Promoterless Expression System	1 kit	VPK-411-DJ
AAV-DJ Helper Free shRNA Expression System (Puro)	1 kit	VPK-412-DJ
AAV-DJ/Helper Free shRNA Expression System (GFP)	1 kit	VPK-413-DJ
scAAV-DJ Helper Free Expression System	1 kit	VPK-430-DJ

AAV-DJ/8 Helper Free Complete Expression Systems

Product Name	Size	Catalog Number
AAV-DJ/8 Helper Free Expression System	1 kit	VPK-410-DJ-8
AAV-DJ/8 Helper Free Bicistronic Expression System (Puro)	1 kit	VPK-415-DJ-8
AAV-DJ/8 Helper Free Bicistronic Expression System (Neo)	1 kit	VPK-416-DJ-8
AAV-DJ/8 Helper Free Bicistronic Expression System (Hygro)	1 kit	VPK-417-DJ-8
AAV-DJ/8 Helper Free Bicistronic Expression System (GFP)	1 kit	VPK-418-DJ-8
AAV-DJ/8 Helper Free Bicistronic Expression System (Blasticidin)	1 kit	VPK-419-DJ-8
AAV-DJ/8 Helper Free Promoterless Expression System	1 kit	VPK-411-DJ-8
AAV-DJ/8 Helper Free shRNA Expression System (Puro)	1 kit	VPK-412-DJ-8
AAV-DJ/8 Helper Free shRNA Expression System (GFP)	1 kit	VPK-413-DJ-8
scAAV-DJ/8 Helper Free Expression System	1 kit	VPK-430-DJ-8

AAV-1 Helper Free Complete Expression Systems

Product Name	Size	Catalog Number
AAV-1 Helper Free Expression System	1 kit	VPK-410-SER1
AAV-1 Helper Free Bicistronic Expression System (Puro)	1 kit	VPK-415-SER1
AAV-1 Helper Free Bicistronic Expression System (Neo)	1 kit	VPK-416-SER1
AAV-1 Helper Free Bicistronic Expression System (Hygro)	1 kit	VPK-417-SER1
AAV-1 Helper Free Bicistronic Expression System (GFP)	1 kit	VPK-418-SER1
AAV-1 Helper Free Bicistronic Expression System (Blasticidin)	1 kit	VPK-419-SER1
AAV-1 Helper Free Promoterless Expression System	1 kit	VPK-411-SER1
AAV-1 Helper Free shRNA Expression System (Puro)	1 kit	VPK-412-SER1
AAV-1 Helper Free shRNA Expression System (GFP)	1 kit	VPK-413-SER1
scAAV-1 Helper Free Expression System	1 kit	VPK-430-SER1

AAV-2 Helper Free Complete Expression Systems

Product Name	Size	Catalog Number
AAV-2 Helper Free Expression System	1 kit	VPK-410-SER2
AAV-2 Helper Free Bicistronic Expression System (Puro)	1 kit	VPK-415-SER2
AAV-2 Helper Free Bicistronic Expression System (Neo)	1 kit	VPK-416-SER2
AAV-2 Helper Free Bicistronic Expression System (Hygro)	1 kit	VPK-417-SER2
AAV-2 Helper Free Bicistronic Expression System (GFP)	1 kit	VPK-418-SER2
AAV-2 Helper Free Bicistronic Expression System (Blasticidin)	1 kit	VPK-419-SER2
AAV-2 Helper Free Promoterless Expression System	1 kit	VPK-411-SER2
AAV-2 Helper Free shRNA Expression System (Puro)	1 kit	VPK-412-SER2
AAV-2 Helper Free shRNA Expression System (GFP)	1 kit	VPK-413-SER2
scAAV-2 Helper Free Expression System	1 kit	VPK-430-SER2

AAV-3 Helper Free Complete Expression Systems

Product Name	Size	Catalog Number
AAV-3 Helper Free Expression System	1 kit	VPK-410-SER3
AAV-3 Helper Free Bicistronic Expression System (Puro)	1 kit	VPK-415-SER3
AAV-3 Helper Free Bicistronic Expression System (Neo)	1 kit	VPK-416-SER3
AAV-3 Helper Free Bicistronic Expression System (Hygro)	1 kit	VPK-417-SER3
AAV-3 Helper Free Bicistronic Expression System (GFP)	1 kit	VPK-418-SER3
AAV-3 Helper Free Bicistronic Expression System (Blasticidin)	1 kit	VPK-419-SER3
AAV-3 Helper Free Promoterless Expression System	1 kit	VPK-411-SER3
AAV-3 Helper Free shRNA Expression System (Puro)	1 kit	VPK-412-SER3
AAV-3 Helper Free shRNA Expression System (GFP)	1 kit	VPK-413-SER3
scAAV-3 Helper Free Expression System	1 kit	VPK-430-SER3

AAV-4 Helper Free Complete Expression Systems

Product Name	Size	Catalog Number
AAV-4 Helper Free Expression System	1 kit	VPK-410-SER4
AAV-4 Helper Free Bicistronic Expression System (Puro)	1 kit	VPK-415-SER4
AAV-4 Helper Free Bicistronic Expression System (Neo)	1 kit	VPK-416-SER4
AAV-4 Helper Free Bicistronic Expression System (Hygro)	1 kit	VPK-417-SER4
AAV-4 Helper Free Bicistronic Expression System (GFP)	1 kit	VPK-418-SER4
AAV-4 Helper Free Bicistronic Expression System (Blasticidin)	1 kit	VPK-419-SER4
AAV-4 Helper Free Promoterless Expression System	1 kit	VPK-411-SER4
AAV-4 Helper Free shRNA Expression System (Puro)	1 kit	VPK-412-SER4
AAV-4 Helper Free shRNA Expression System (GFP)	1 kit	VPK-413-SER4
scAAV-4 Helper Free Expression System	1 kit	VPK-430-SER4

AAV-5 Helper Free Complete Expression Systems

Product Name	Size	Catalog Number
AAV-5 Helper Free Expression System	1 kit	VPK-410-SER5
AAV-5 Helper Free Bicistronic Expression System (Puro)	1 kit	VPK-415-SER5
AAV-5 Helper Free Bicistronic Expression System (Neo)	1 kit	VPK-416-SER5
AAV-5 Helper Free Bicistronic Expression System (Hygro)	1 kit	VPK-417-SER5
AAV-5 Helper Free Bicistronic Expression System (GFP)	1 kit	VPK-418-SER5
AAV-5 Helper Free Bicistronic Expression System (Blasticidin)	1 kit	VPK-419-SER5
AAV-5 Helper Free Promoterless Expression System	1 kit	VPK-411-SER5
AAV-5 Helper Free shRNA Expression System (Puro)	1 kit	VPK-412-SER5
AAV-5 Helper Free shRNA Expression System (GFP)	1 kit	VPK-413-SER5
scAAV-5 Helper Free Expression System	1 kit	VPK-430-SER5

AAV-6 Helper Free Complete Expression Systems

Product Name	Size	Catalog Number
AAV-6 Helper Free Expression System	1 kit	VPK-410-SER6
AAV-6 Helper Free Bicistronic Expression System (Puro)	1 kit	VPK-415-SER6
AAV-6 Helper Free Bicistronic Expression System (Neo)	1 kit	VPK-416-SER6
AAV-6 Helper Free Bicistronic Expression System (Hygro)	1 kit	VPK-417-SER6
AAV-6 Helper Free Bicistronic Expression System (GFP)	1 kit	VPK-418-SER6
AAV-6 Helper Free Bicistronic Expression System (Blasticidin)	1 kit	VPK-419-SER6
AAV-6 Helper Free Promoterless Expression System	1 kit	VPK-411-SER6
AAV-6 Helper Free shRNA Expression System (Puro)	1 kit	VPK-412-SER6
AAV-6 Helper Free shRNA Expression System (GFP)	1 kit	VPK-413-SER6
scAAV-6 Helper Free Expression System	1 kit	VPK-430-SER6

AAV Helper Free Packaging Systems

AAV Helper Free Packaging Systems contain everything found in the Complete Expression Systems, with the exception of the AAV expression vector. This is an ideal choice if you already have an AAV construct containing your gene of interest.

The following vectors are included with each AAV Helper Free Packaging System:

- pHelper Plasmid
- Rep-Cap Plasmid (serotype specific)
- GFP Control Vector

All plasmids are provided individually, not as a packaging mixture.

Recent Product Citations

1. Iwata, R. et al. (2015). Developmental RacGAP α 2-chimaerin signaling is a determinant of the morphological features of dendritic spines in adulthood. *J. Neurosci.* **35**:13728-13744. (VPK-400-DJ-8)
2. Friedland, A.E. et al. (2015). Characterization of *Staphylococcus aureus* Cas9: a smaller Cas9 for all-in-one adeno-associated virus delivery and paired nickase applications. *Genome Biol.* **16**:257. (VPK-402)

- If you need an AAV expression vector as well as an AAV packaging system, choose one of our **AAV Complete Expression Systems: p. 42-45**
- If you want to make AAV from more than one serotype but don't want to order multiple packaging systems, choose one AAV Packaging System from this list and then add additional **AAV Rep-Cap Plasmids: p. 46**

Product Name	Size	Catalog Number
AAV-DJ Helper Free Packaging System	1 kit	VPK-400-DJ
AAV-DJ/8 Helper Free Packaging System	1 kit	VPK-400-DJ-8
AAV-1 Helper Free Packaging System	1 kit	VPK-401
AAV-2 Helper Free Packaging System	1 kit	VPK-402
AAV-3 Helper Free Packaging System	1 kit	VPK-403
AAV-4 Helper Free Packaging System	1 kit	VPK-404
AAV-5 Helper Free Packaging System	1 kit	VPK-405
AAV-6 Helper Free Packaging System	1 kit	VPK-406

AAV Rep-Cap Plasmids (Serotype-Specific)

AAV Rep-Cap plasmids allow you to make recombinant AAV of a specific serotype. These plasmids are ideal if you already have an AAV packaging system for a different serotype. Just substitute one of these plasmids into your AAV Helper Free Packaging System or Expression System.

Recent Product Citations

- Blackburn, J. et al. (2015). Damage-inducible intragenic demethylation of the human TP53 tumor suppressor gene is associated with transcription from an alternative intronic promoter. *Mol. Carcinog.* 10.1002/mc.22441. (VPK-420-DJ)
- Kato, H.E. et al. (2015). Atomistic design of microbial opsin-based blue-shifted optogenetics tools. *Nat. Commun.* 6:7177. (VPK-420-DJ)
- Holehonnur, R. et al. (2014). Adeno-associated viral serotypes produce differing titers and differentially transduce neurons within the rat basal and lateral amygdala. *BMC Neurosci.* 15:28. (VPK-420-DJ, VPK-420-DJ-8)

Product Name	Catalog Number
pAAV-DJ Vector	VPK-420-DJ
pAAV-DJ/8 Vector	VPK-420-DJ-8
pAAV-RC1 Vector	VPK-421
pAAV-RC2 Vector	VPK-422

Product Name	Catalog Number
pAAV-RC3 Vector	VPK-423
pAAV-RC4 Vector	VPK-424
pAAV-RC5 Vector	VPK-425
pAAV-RC6 Vector	VPK-426

AAV Expression Vectors

Each of our AAV Expression Vectors may be used with any of our AAV Helper Free Systems, regardless of AAV serotype. Choose one of these vectors when you already have an AAV Packaging System but may want to use a different promoter or selection marker.

Recent Product Citations

- Tang, F.L. et al. (2015). VPS35 deficiency or mutation causes dopaminergic neuronal loss by impairing mitochondrial fusion and function. *Cell Rep.* 12:1631-1643. (VPK-410, VPK-413)
- Orabi, A.I. et al. (2015). Dynamic imaging of pancreatic NF- κ B activation in live mice using AAV infusion and bioluminescence. *J. Biol. Chem.* 10.1074/jbc.M115.647933. (VPK-410)
- Fyk-Kolodziej, B. et al. (2014). Marking cells with infrared fluorescent proteins to preserve photoresponsiveness in the retina. *Biotechniques* 57:245-253. (VPK-410)
- Lutz, D. et al. (2014). Myelin basic protein cleaves cell adhesion molecule L1 and promotes neuritogenesis and cell survival. *J. Biol. Chem.* 289:13503-13518. (VPK-410)
- Xiang, J. et al. (2015). Postnatal loss of Hap1 reduces hippocampal neurogenesis and causes adult depressive-like behavior in mice. *PLoS Genet.* 11:e1005175. (VPK-411)
- Nishikawa, T. et al. (2015). Resetting the transcription factor network reverses terminal chronic hepatic failure. *J. Clin. Invest.* 10.1172/JCI73137. (VPK-418)
- Huckstepp, R.T. et al. (2015). Role of parafacial nuclei in control of breathing in adult rats. *J. Neurosci.* 35:1052-1067. (VPK-418)
- Krol, J. et al. (2015). A network comprising short and long noncoding RNAs and RNA helicase controls mouse retina architecture. *Nat. Commun.* 6:7305. (VPK-430)

Product Name	Cloning Capacity	Size	Catalog Number
pAAV-MCS Expression Vector	3 kb	10 μ g	VPK-410
pAAV-IRES-Puro Expression Vector	1.8 kb	10 μ g	VPK-415
pAAV-IRES-Neo Expression Vector	1.6 kb	10 μ g	VPK-416
pAAV-IRES-Hygro Expression Vector	1.4 kb	10 μ g	VPK-417
pAAV-IRES-GFP Expression Vector	1.7 kb	10 μ g	VPK-418
pAAV-IRES-Bsd Expression Vector	2 kb	10 μ g	VPK-419
pAAV-MCS Promoterless Expression Vector	3.9 kb	10 μ g	VPK-411
pAAV-U6-Puro Expression Vector	2.2 kb	10 μ g	VPK-412
pAAV-U6-GFP Expression Vector	2.1 kb	10 μ g	VPK-413
pscAAV-MCS Expression Vector	1.5 kb	10 μ g	VPK-430

AAV Control Plasmids

Our AAV control plasmids are useful for making a transduction control virus. Choose one of these constructs to pair with an AAV Packaging System. The vectors are provided as purified DNA.

Product Name	Size	Catalog Number
pAAV-GFP Control Vector	10 µg	AAV-400
pAAV-Cre Control Vector	10 µg	AAV-401
pAAV-LacZ Control Vector	10 µg	AAV-402
pscAAV-GFP Control Vector	10 µg	AAV-410

AAV Premade Control Viruses

Our AAV premade control viruses provide a convenient tool for measuring transduction efficiency into a target cell. Each virus contains a reporter gene, with the exception of the Null Control virus.

All AAV premade viruses are provided at a concentration of 1×10^{12} GC/mL.

Product Name	Size	Catalog Number
AAV1-GFP Control Virus	50 µL	AAV-301
AAV2 Null Control Virus	50 µL	AAV-300
AAV2-Cre Control Virus	50 µL	AAV-310
AAV2-GFP Control Virus	50 µL	AAV-302
AAV2-Luc Control Virus	50 µL	AAV-320
AAV3-GFP Control Virus	50 µL	AAV-303
AAV5-GFP Control Virus	50 µL	AAV-305
AAV6-GFP Control Virus	50 µL	AAV-306

293AAV Cell Line

Our 293AAV cell line is selected from the parental 293 cell line for larger surface area, flattened morphology, and firmer attachment to culture plates, resulting in production of higher yields of AAV.

Recent Product Citation

Lutz, D. et al. (2014). Myelin basic protein cleaves cell adhesion molecule L1 and promotes neuritogenesis and cell survival. *J. Biol. Chem.* **289**:13503-13518.

Product Name	Size	Catalog Number
293AAV Cell Line	1×10^6 cells	AAV-100


ViraBind™ AAV Purification Kits

Purification of AAV via ultracentrifugation can be tedious and time-consuming, and may result in low yields. ViraBind™ AAV Purification Kits use a one-step proprietary matrix followed by further purification and concentration using a centrifugal concentrator. The result is a higher AAV yield with high purity in a fraction of the time. Kits are available in two sizes:

- The Standard kit can purify up to two 10cm dishes per prep
- The Mega kit can purify up to ten 15cm dishes per prep

ViraBind™ Kits are suitable for AAV-2 or AAV-DJ; they are not compatible with other AAV serotypes.


- **High Purity:** No contamination bands as seen on SDS gel
- **Fast Results:** Obtain purified virus in about 3 hours
- **High Yields:** Recovery rate >60%



Purification Procedure for the ViraBind™ AAV Purification Kit.

Recent Product Citations

1. Tsuneoka, Y. et al. (2015). Distinct preoptic-BST nuclei dissociate paternal and infanticidal behavior in mice. *EMBO J.* 10.15252/embj.201591942. (VPK-140)
2. Nishikawa, T. et al. (2015). Resetting the transcription factor network reverses terminal chronic hepatic failure. *J. Clin. Invest.* 10.1172/JCI73137. (VPK-140)
3. Stankowska, D.L. et al. (2015). Neuroprotective effects of transcription factor Brn3b in an ocular hypertension rat model of glaucoma. *Invest. Ophthalmol. Vis. Sci.* 56:893-907. (VPK-140)
4. Moshiri, F. et al. (2014). Inhibiting the oncogenic mir-221 by microRNA sponge: toward microRNA-based therapeutics for hepatocellular carcinoma. *Gastroenterol. Hepatol. Bed Bench* 7:43-54. (VPK-140)
5. Inaba, Y. et al. (2014). Gadd34 regulates liver regeneration in hepatic steatosis. *Hepatology* 6:1343-1356. (VPK-140)
6. Rodriguez, J.P. et al. (2014). Abrogation of β -catenin signaling in oligodendrocyte precursor cells reduces glial scarring and promotes axon in the medial prefrontal cortex of rodents. *J. Neurosci.* 30:15007-15018. (VPK-140)
7. Uchida, S. et al. (2010). Early life stress enhances behavioral vulnerability to stress through the activation of REST4-mediated gene transcription in the medial prefrontal cortex of rodents. *J. Neurosci.* 30:15007-15018. (VPK-140)
8. Nazari, M. et al. (2014). AAV2-mediated follistatin overexpression induces ovine primary myoblasts proliferation. *BMC Biotechnol.* 14:87. (VPK-141)




Electrophoretic Profile of Purified AAV2-GFP.

Product Name	Capacity/Prep	Size	Catalog Number
ViraBind™ AAV Purification Kit	Two 10-cm dishes	10 Preps	VPK-140
ViraBind™ AAV Purification Mega Kit	Ten 15-cm dishes	2 Preps	VPK-141
		10 Preps	VPK-141-5

QuickTiter™ AAV Quantitation Kit

Traditional AAV Quantitation by dot blot can be tedious, time consuming, and suffer from high inter-assay variability. Our QuickTiter™ AAV Quantitation Kit uses a proprietary technology to quantify AAV nucleic acid content of unpurified AAV-2 or AAV-DJ, or from purified AAV of any serotype.



AAV-2 DNA Standard Curve. The QuickTiter™ AAV-2 DNA Standard was diluted as described in the assay protocol. Fluorescence was measured on a SpectraMax Gemini XS Fluorometer (Molecular Devices) with a 485/538 nm filter set and 530 nm cutoff.

- **Fast Results:** Obtain purified virus in less than 2 hours
- **High Sensitivity:** Limit of detection 1×10^9 GC/mL from unpurified supernatant or 5×10^{10} GC/mL from purified AAV

Recent Product Citations

1. Orabi, A.I. et al. (2015). Dynamic imaging of pancreatic NF- κ B activation in live mice using AAV infusion and bioluminescence. *J. Biol. Chem.* 10.1074/jbc.M115.647933.
2. Oh, S.M. et al. (2015). Combined Nurr1 and Foxa2 roles in the therapy of Parkinson's disease. *EMBO Mol. Med.* 10.15252/emmm.201404610.
3. Stankowska, D.L. et al. (2015). Neuroprotective effects of transcription factor Brn3b in an ocular hypertension rat model of glaucoma. *Invest. Ophthalmol. Vis. Sci.* 56:893-907.
4. Li, Y. et al. (2014). Assembly and validation of versatile transcription activator-like effector libraries. *Sci. Rep.* 4:4857.
5. Paydar, A. et al. (2014). Extrasynaptic GABAA receptors in mediodorsal thalamic nucleus modulate fear extinction learning. *Mol. Brain* 7:39.
6. Vinnikov, I.A. et al. (2014). Hypothalamic miR-103 protects from hyperphagic obesity in mice. *J. Neurosci.* 34:10659-10674.
7. Ma, J. et al. (2013). RNA interference-mediated silencing of Atp6i prevents both periapical bone erosion and inflammation in the mouse model of endodontic disease. *Infect. Immun.* 81:1021-1030.
8. Tao, P. et al. (2013). In vitro and in vivo delivery of genes and proteins using the bacteriophage T4 DNA packaging machine. *PNAS* 10.1073/pnas.1300867110.

Product Name	Capacity/Prep	Size	Catalog Number
QuickTiter™ AAV Quantitation Kit	Fluorometric	20 Assays	VPK-145

ViraDuctin™ AAV Transduction Reagent

Successful gene expression studies using AAV depend on high transduction efficiencies into host cells. Infection rates appear to be highest in S-phase cells, which can account for a very small fraction of a cell population.

Our ViraDuctin™ AAV Transduction Reagent can significantly increase the transduction efficiency of AAV vectors in both dividing and non-dividing cells. Increases are greatest in non-dividing cells, but even cells in S-phase show a noticeable increase in transduction efficiencies.

- **Higher Efficiencies:** Significantly increase rate of infection of host cells
- **Low Toxicity:** No noticeable effect on cell viability
- **Universal:** Suitable for use with both dividing and non-dividing cells

Recent Product Citation

Nazari, M. et al. (2014). AAV2-mediated follistatin overexpression induces ovine primary myoblasts proliferation. *BMC Biotechnol.* 14:87.

Product Name	Size*	Catalog Number
ViraDuctin™ AAV Transduction Reagent	10 Transductions	AAV-200
	50 Transductions	AAV-201

*Number of transductions performed in 35mm culture dishes. May be modified for use in culture plates or larger dishes. See product insert.

Adenoviral Expression Kits & Reagents


Recombinant adenoviruses are excellent tools for introducing genetic material into host cells, since they can infect a variety of mammalian cell types with high efficiency. They remain epichromosomal upon infection, so they are only suitable for transient gene expression. We offer a complete workflow solution to your adenoviral expression studies:

- Premade Recombinant Adenoviruses
- Viral Expression Systems
- Viral Packaging Cell Line
- Purification Kits
- Quantitation / Titer Kits
- Transduction Reagents

Premade Recombinant Adenoviruses

Don't have time to make your own adenovirus? Are you studying the expression of multiple genes? Rely on our premade recombinant adenoviruses that already contain a gene of interest. All of Cell Biolabs' premade recombinant adenoviruses are provided as 50 µl aliquots at a concentration of 1×10^{11} viral particles/mL in TBS with 10% glycerol.

Angiogenesis



Blood Vessel Formation After 3 Days. Purified Ad-Null or Ad-VEGF viruses were applied to a 10-day old CAM (chick chorioallantoic membrane). Results were visualized by stereomicroscope.

Recent Product Citations

1. Jiang, Y.Z. et al. (2014). Distinct roles of HIF1A in endothelial adaptations to physiological and ambient oxygen. *Mol. Cell Endocrinol.* **391**:60-67. (ADV-100)
2. Kelber, J.A. et al. (2012). KRas induces a Src/PEAK1/ErbB2 kinase amplification loop that drives metastatic growth and therapy resistance in pancreatic cancer. *Cancer Res.* **72**:2554-2564. (ADV-101)
3. Qiu, X. et al. (2012). Combined strategy of mesenchymal stem cell injection with vascular endothelial growth factor gene therapy for the treatment of diabetes-associated erectile dysfunction. *J. Androl.* **33**:37-44. (ADV-101)
4. Stoletov, K. et al. (2010). Visualizing extravasation dynamics of metastatic tumor cells. *J. Cell Sci.* **123**:2332-2341. (ADV-101)

Target Name	Catalog Number
HIF-1 α	ADV-100
VEGF	ADV-101

Controls and Reporter Genes

Recent Product Citations

1. Wang, Y. et al. (2015). MiR-31 downregulation protects against cardiac ischemia/reperfusion injury by targeting protein kinase C epsilon (PKC ϵ) directly. *Cell Physiol. Biochem.* **36**:179-190. (ADV-001)
2. Haidari, M. et al. (2014). Disruption of endothelial adherens junctions by high glucose is mediated by protein kinase C- β -dependent vascular endothelial cadherin tyrosine phosphorylation. *Circulation.* **130**:112. (ADV-001 and ADV-004)
3. Kaneshiro, S. et al (2015). MEK5 suppresses osteoblastic differentiation. *Biochem. Biophys. Res. Commun.* **463**:1016-1019. (ADV-002)
4. Choi, J.M. et al (2015). HepG2 cells as an in vitro model for evaluation of cytochrome P450 induction by xenobiotics. *Arch. Pharm. Res.* **38**:691-704. (ADV-002 and ADV-005)
5. Lampert, F.M. et al. (2015). Overexpression of Hif-1a in mesenchymal stem cells affects cell-autonomous angiogenic and osteogenic parameters. *J. Cell Biochem.* **120**:2536-2546. (ADV-004)
6. Jiang, Y.Z. et al. (2014). Distinct roles of HIF1A in endothelial adaptations to physiological and ambient oxygen. *Mol. Cell Endocrinol.* **391**:60-67. (ADV-004)
7. Salvati, E. et al. (2014). Evidence for G-quadruplex in the promoter of VEGFR-2 and its targeting to inhibit tumor angiogenesis. *Nucleic Acids Res.* **42**:2945-2957. (ADV-004)
8. Wu, Y. et al. (2012). ERK5 regulates glucose-induced increased fibronectin production in the endothelial cells and in the retina in diabetes. *Invest. Ophthalmol. Vis. Sci.* **53**:8405-8413. (ADV-004)
9. Choi, S. et al. (2014). MMP9 processing of HSPB1 regulates tumor progression. *PLoS One* **9**:e85509. (ADV-005)
10. Kato, H. et al. (2011). Wnt/ β -Catenin pathway in podocytes integrates cell adhesion, differentiation and survival. *J. Biol. Chem.* **286**:26003-26015. (ADV-005)

Target Name	Catalog Number
Null Control (No gene)	ADV-001
β -Galactosidase	ADV-002
Cre	ADV-005
Firefly Luciferase	ADV-008
GFP	ADV-004
SEAP (Secretory Alkaline Phosphatase)	ADV-003

Premade Recombinant Adenoviruses, continued

MAP Kinase Signaling

Recent Product Citations

1. Harbrecht, B.G. et al. (2012). Insulin inhibits hepatocyte iNOS expression induced by cytokines by an Akt-dependent mechanism. *Am. J. Physiol. Gastrointest. Liver Physiol* **302**:G116-G122. (ADV-105)
2. Monick, M. et al. (2008). Constitutive ERK MAPK activity regulates macrophage ATP production and mitochondrial integrity. *J. Immunol.* **180**:7485-7496. (ADV-112, ADV-113, ADV-118, ADV-119)
3. Li, M.Y. et al. (2014). Curcumin inhibits 19-kDa lipoprotein of *Mycobacterium tuberculosis* induced macrophage apoptosis via regulation of the JNK pathway. *Biochem. Biophys. Res. Commun.* **446**:626-632. (ADV-115)
4. Jiang, S. et al. (2014). Regulation of hepatic insulin receptor activity following injury. *Am. J. Physiol. Gastrointest. Liver Physiol.* **306**:G886-G892. (ADV-115)
5. Jiang, S. et al. (2011). Role of inhibitory kB kinase and c-Jun NH2-terminal kinase in the development of hepatic insulin resistance in critical illness diabetes. *Am. J. Physiol. Gastrointest. Liver Physiol.* **301**:G454-G463. (ADV-115)
6. Zhao, P. et al. (2015). Filamin A (FLNA) modulates chemosensitivity to docetaxel in triple-negative breast cancer through the MAPK/ERK pathway. *Tumor Biol.* 10.1007/s13277-015-4357-3. (ADV-118)
7. Aissaoui, H. et al. (2015). MDA-9/syntenin is essential for factor VIIa-induced signaling, migration, and metastasis in melanoma cells. *J. Biol. Chem.* **290**:3333-3348. (ADV-118)
8. Zhang, Z. et al. (2013). MEK inhibition leads to lysosome-mediated Na⁺/I⁻ symporter protein degradation in human breast cancer cells. *Endocr. Relat. Cancer* **20**:241-250. (ADV-118)
9. Goc, A. et al. (2015). p70 S6-kinase mediates the cooperation between Akt1 and MEK1 pathways in fibroblast-mediated extracellular matrix remodeling. *Biochim. Biophys. Acta.* **1853**:1626-1635. (ADV-118 and ADV-119)
10. Matsushita, T. et al. (2009). FGFR3 promotes synchondrosis closure and fusion of ossification centers through the MAPK pathway. *Hum. Mol. Genet.* **18**:227-240. (ADV-119)
11. Tan, S.H. et al. (2009). Regulation of cell proliferation and migration by TAK1 via transcriptional control of von Hippel-Lindau tumor suppressor. *J. Biol. Chem.* **284**:18047-18058. (ADV-128)
12. Kaneshiro, S. et al (2015). MEK5 suppresses osteoblastic differentiation. *Biochem. Biophys. Res. Commun.* 10.1016/j.bbrc.2015.05.035. (ADV-129)
13. Stankiewicz, T.R. et al. (2015). Neuronal apoptosis induced by selective inhibition of Rac GTPase versus global suppression of Rho family GTPases is mediated by alterations in distinct mitogen-activated protein kinase signaling cascades. *J. Biol. Chem.* **290**:9363-9376. (ADV-129 and ADV-131)
14. Wu, Y. et al. (2012). ERK5 regulates glucose-induced increased fibronectin production in the endothelial cells and in the retina in diabetes. *Invest. Ophthalmol. Vis. Sci.* **53**:8405-8413. (ADV-130)
15. Zuo, Y. et al. (2015). Modulation fo ERK5 is a novel mechanism by which Cdc42 regulates migration of breast cancer cells. *J. Cell Biochem.* **116**:124-132. (ADV-131)
16. Ozcan, L. et al. (2015). Treatment of obese insulin-resistant mice with an allosteric MAPKAPK2/3 inhibitor lowers blood glucose and improves insulin sensitivity. *Diabetes* **64**:3396-3405. (ADV-138)
17. Taniguchi, C. et al. (2007). The p85a regulatory subunit of phosphoinositide 3-kinase potentiates c-Jun N-terminal kinase-mediated insulin resistance. *Mol. Cell Biol.* **27**:2830-2840. (ADV-161)

Target Name	Catalog Number
ERK2	ADV-112
ERK2 (Dominant Negative)	ADV-113
Interferon- γ	ADV-103
JNK1	ADV-114
JNK1 (Dominant Negative)	ADV-115
MAPKAPK2	ADV-137
MAPKAPK2 (Dominant Negative)	ADV-138
MAPKAPK2 (Constitutively Active)	ADV-139
MEK1 (Dominant Negative)	ADV-118
MEK1 (Constitutively Active)	ADV-119
MEK5	ADV-129
MEK5 (Dominant Negative)	ADV-130
MEK5 (Constitutively Active)	ADV-131
MEKK1	ADV-135
MEKK1 (Dominant Negative)	ADV-136
MKK3	ADV-120
MKK3 (Dominant Negative)	ADV-121
MKK3 (Constitutively Active)	ADV-122
MKK4 (Dominant Negative)	ADV-160
MKK4 (Constitutively Active)	ADV-161
MKK6	ADV-123
MKK6 (Dominant Negative)	ADV-124
MKK6 (Constitutively Active)	ADV-125
MKK7	ADV-126
MKK7 (Dominant Negative)	ADV-127
MKK7 (Constitutively Active)	ADV-128
p38 α	ADV-104
p38 α (Dominant Negative)	ADV-105
Raf1	ADV-132
Raf1 (Dominant Negative)	ADV-133
Raf1 (Constitutively Active)	ADV-134

Premade Recombinant Adenoviruses, continued

Cytoskeleton Regulation / Small GTPase

Recent Product Citations

1. Aissaoui, H. et al. (2015). MDA-9/syntenin is essential for factor VIIa-induced signaling, migration, and metastasis in melanoma cells. *J. Biol. Chem.* **290**:3333-3348. (ADV-149, ADV-152)
2. Mao, Y. et al. (2012). Essential diurnal Rac1 activation during retinal phagocytosis requires $\alpha\beta 5$ integrin but not tyrosine kinases focal adhesion kinase or Mer tyrosine kinase. *Mol. Cell Biol.* **23**:1104-1114. (ADV-150)
3. Thomas, M.A. et al. (2009). E4orf1 limits the oncolytic potential of the E1B-55K-deleted adenovirus. *J. Virol.* **83**:2406-2416. (ADV-150)
4. Yu, W.-M. et al. (2009). Laminin is required for Schwann cell morphogenesis. *J. Cell Sci.* **122**:929-936. (ADV-150, ADV-153, ADV-154)
5. Salvati, E. et al. (2014). Evidence for G-quadruplex in the promoter of VEGFR-2 and its targeting to inhibit tumor angiogenesis. *Nucleic Acids Res.* **42**:2945-2957. (ADV-151, ADV-157)
6. Cheng, Z.-J. et al. (2010). Co-regulation of caveolar and Cdc42-dependent fluid phase endocytosis by phosphocaveolin-1. *J. Biol. Chem.* **285**:15119-15125. (ADV-153)
7. Neal M. et al. (2013). A critical role for TLR4 induction of autophagy in the regulation of enterocyte migration and the pathogenesis of necrotizing enterocolitis. *J. Immunol.* **190**:3541-3551. (ADV-156, ADV-157)

Target Name	Catalog Number
Cdc42	ADV-152
Cdc42 L61 (Constitutively Active)	ADV-154
Cdc42 N17 (Dominant Negative)	ADV-153
Rac1	ADV-149
Rac1 L61 (Constitutively Active)	ADV-151
Rac1 N17 (Dominant Negative)	ADV-150
Ras N17 (Dominant Negative)	ADV-145
Ras V12 (Constitutively Active)	ADV-146
Ras V12C40	ADV-148
Ras V12S35	ADV-147
Rho L63 (Constitutively Active)	ADV-157
Rho N19 (Dominant Negative)	ADV-156
SDF-1 α	ADV-210

Cell Cycle & Transcription Regulation

Recent Product Citation

Nguyen, N. et al. (2014). Mitsugumin 53 (MG53) ligase ubiquitinates focal adhesion kinase during skeletal myogenesis. *J. Biol. Chem.* **289**:3209-3216. (ADV-508)

Target Name	Catalog Number
MyoD	ADV-508
p53	ADV-501
p53 (Temperature Sensitive Mutant)	ADV-502

NF κ B Signaling**Recent Product Citation**

Wang, Y. et al. (2015). MiR-31 downregulation protects against cardiac ischemia/reperfusion injury by targeting protein kinase C epsilon (PKC ϵ) directly. *Cell Physiol. Biochem.* **36**:179-190. (ADV-302)

Target Name	Catalog Number
I κ B- α S32A (Dominant Negative)	ADV-302
IKK- β	ADV-305
IKK- β (Dominant Negative)	ADV-303

293AD Cell Line for Adenoviral Packaging and Amplification

The 293AD cell line is ideal for packaging and amplifying adenovirus. The cell line is derived from the parental 293 cell line, but has been specifically selected for adenovirus applications and offers advantages over conventional 293 cells: flattened morphology, firm attachment to culture plates, and a larger surface area for superior transfection and greater viral yields.

Recent Product Citations

1. Bae, E.J. et al. (2015). Cell models to study cell-to-cell transmission of α -synuclein. *Methods Mol. Biol.* **1345**:291-298.
2. Sugiyama, K. et al. (2014). Expression of the miR200 family of microRNAs in mesothelial cells suppresses the dissemination of ovarian cancer cells. *Mol. Cancer Ther.* **13**:2081-2091.
3. Peng, D. et al. (2014). Glutathione peroxidase 7 has potential tumour suppressor functions that are silenced by location-specific methylation in oesophageal adenocarcinoma. *Gut* **63**:540-551.

Product Name	Size	Catalog Number
293AD Cell Line	1 x 10 ⁶ Cells	AD-100


RAPAd® Adenoviral Expression Systems

Compared to other adenoviral expression systems, RAPAd® Adenoviral Expression Systems produce recombinant adenovirus in a much shorter time (about 2-3 weeks) with a substantial reduction in wild-type adenovirus. The RAPAd systems use a backbone vector from which the 5' ITR, packaging signal and E1 sequences have been removed. Additionally, serial amplification of the recombinant adenovirus does not increase the level of replication-competent adenovirus.

Recent Product Citations

- Li, P. et al. (2013). MicroRNA-663 regulates human vascular smooth muscle cell phenotypic switch and vascular neointimal formation. *Circ. Res.* **113**:1117-1127. (VPK-250)
- Bae, E.J. et al. (2015). Cell models to study cell-to-cell transmission of α -synuclein. *Methods Mol. Biol.* **1345**:291-298. (VPK-252)
- Suchanek, A.L. and Salati, L.M. (2015). Construction and evaluation of an adenoviral vector for the liver-specific expression of the serine/arginine rich splicing factor, SRSF3. *Plasmid* 10.1016/j.plasmid.2015.07.004. (VPK-252)
- Li, M. et al. (2014). Bisphenol AF-induced endogenous transcription is mediated by ER α and ERK1/2 activation in human breast cancer cells. *PLoS One* **9**:e94725. (VPK-252)
- Kim, S.C. et al. (2014). All-trans-retinoic acid ameliorates hepatic steatosis in mice by a novel transcriptional cascade. *Hepatology* **59**:1750-1760. (VPK-252)
- Mohan, R. et al. (2015). Differentially expressed microRNA-483 confers distinct functions in pancreatic beta- and alpha-cells. *J. Biol. Chem.* 10.1074/jbc.M115.650705. (VPK-253)
- Sugiyama, K. et al. (2014). Expression of the miR200 family of microRNAs in mesothelial cells suppresses the dissemination of ovarian cancer cells. *Mol. Cancer Ther.* **13**:2081-2091. (VPK-253)
- Feng, J. et al. (2015). SIRT6 suppresses glioma cell growth via induction of apoptosis, inhibition of oxidative stress and suppression of JAK2/STAT3 signaling pathway activation. *Oncol Rep.* 10.3892/or.2015.4477. (VPK-254)
- Lozic, M. et al. (2015). Overexpression of oxytocin receptors in the hypothalamic PVN increases baroreceptor reflex sensitivity and buffers BP variability in conscious rats. *Br. J. Pharmacol.* **171**:4385-4398. (VPK-254)
- Brunton, P.J. et al. (2015). 5 α -reduced neurosteroids sex-dependently reverse central prenatal programming of neuroendocrine stress responses in rats. *J. Neurosci.* **35**:666-677. (VPK-254)

- **Virtually No Wild-Type Virus:** Backbone vector engineered to produce <300 wild-type plaques per 10⁹ particles, compared with 10⁴-10⁶ WT plaques per 10⁹ VP with most other methods
- **Faster Production:** Virus generated in 2-3 weeks compared to a few months with traditional methods
- **7 Complete Systems:** Choose CMV or RSV for gene expression, EF-1 for miRNA expression, U6 for shRNA, or clone your own promoter along with your gene of interest using our Universal system




Adenovirus Production using the RAPAd® Adenoviral Expression System.

Product Name	Promoter	Size	Catalog Number
RAPAd® Universal Adenoviral Expression System	None	1 Kit	VPK-250
RAPAd® RSV Adenoviral Expression System	RSV	1 Kit	VPK-251
RAPAd® CMV Adenoviral Expression System	CMV	1 Kit	VPK-252
RAPAd® Bicistronic Adenoviral Expression System (GFP)	CMV	1 Kit	VPK-254
RAPAd® miRNA Adenoviral Expression System	EF-1	1 Kit	VPK-253
RAPAd® shRNA Adenoviral Expression System (Puro)	U6	1 Kit	VPK-255
RAPAd® shRNA Adenoviral Expression System (GFP)	U6	1 Kit	VPK-256

ViraBind™ Adenovirus Purification Kits

Purification of viruses via cesium chloride (CsCl) ultracentrifugation procedures can be tedious and time-consuming. ViraBind™ Adenovirus Purification Kits use an efficient system for quick adenoviral purification with high recovery. No ultracentrifugation is required. Kits use either a spin column or syringe filter for high purity adenovirus (see selection guide).

- **High Viral Yield:** >90% recovery
- **High Quality:** Provides quality of CsCl procedures, but in much less time
- **Faster Results:** 30 minutes (1-2 hrs for Mega kit)
- **User-Friendly Protocol:** No gradient preparation or ultracentrifugation steps



Purification of Recombinant Ad-β-Gal with the ViraBind™ Adenovirus Purification Kit. Ad-β-Gal was purified according to the assay protocol. Each purification fraction was used to infect A549 cells in a 12-well plate. After 48 hr, cells were scored using our β-Galactosidase Staining Kit (p. 114).

Recent Product Citations

1. Morris, S.J. et al. (2015). Laboratory-scale production of replication-deficient adenovirus vectored vaccines. *Vaccine Technol. Vet Viral Disease* 10.1007/978-1-4939-3008-1_8. (VPK-099)
2. Muller, J. et al. (2015). TROM21, a negative modulator of LFG in breast carcinoma MDA-MB-231 cells in vitro. *Int. J. Oncol.* 47:1634-1646. (VPK-099)
3. Zhu, L. et al. (2015). Inhibition of porcine reproductive and respiratory syndrome virus replication with exosome-transferred artificial microRNA targeting the 3' untranslated region. *J. Virol. Methods* 223:61-68. (VPK-099)
4. Boehme, P. et al. (2015). Standard free droplet digital polymerase chain reaction as a new tool for the quality control of high-capacity adenoviral vectors in small-scale preparations. *Hum. Gene Ther. Methods* 26:25-34. (VPK-099)
5. Kumar, A. et al. (2015). Immune responses against hepatitis C virus genotype 3a virus-like particles in mice: a novel VLP prime-adenovirus boost strategy. *Vaccine* 10.1016/j.vaccine.2015.11.061. (VPK-100)
6. Gibot, L. et al. (2015). Cell-based approach for 3D reconstruction of lymphatic capillaries in vitro reveals distinct functions of HGF and VEGF-C in lymphangiogenesis. *Biomaterials* 10.1016/j.biomaterials.2015.11.027. (VPK-100)
7. Betz, V.M. et al. (2015). Gene-activated fat grafts for the repair of spinal cord injury: a pilot study. *Acta Neurochir.* 10.1007/s00701-0154-2626-y. (VPK-100)
8. Garcia-Pascual, C.M. et al. (2015). Evaluation of the potential therapeutic effects of a double-stranded RNA mimic complexed with polycations in an experimental mouse model of endometriosis. *Fertil. Steril.* 10.1016/j.fertnstert.2015.07.1147. (VPK-100)
9. Tu, X. et al. (2015). microRNA-30 protect against CCl4-induced liver fibrosis by attenuating TGF-β signaling in hepatic stellate cells. *Toxicol. Sci.* 10.1093/toxsci/kfv081. (VPK-100)
10. Miklavc, P. et al. (2015). Actin depolymerisation and crosslinking joint forces with myosin II to contract actin coats on fused secretory vesicles. *J. Cell Sci.* 128:1193-1203. (VPK-100)
11. Westermeier, F. et al. (2015). Insulin requires normal expression and signaling of insulin receptor A to a reverse gestational diabetes-reduced adenosine transport in human umbilical vein endothelium. *FASEB J.* 29:37-49. (VPK-100)

Selection Guide for ViraBind™ Adenovirus Purification Kits


	ViraBind™ Adenovirus Miniprep Kit	ViraBind™ Adenovirus Purification Kit
Purification Method	Spin column	Syringe filter
Purification Time	30 minutes	30 minutes
Capacity/Prep (Viral Particles)	1 x 10 ¹¹ VP	2.5 x 10 ¹² VP
Capacity/Prep (Supernatant)	One T75 flask or one 10cm dish	Four T75 flasks

Product Name	Capacity/Prep	Size	Catalog Number
ViraBind™ Adenovirus Miniprep Kit	1 x 10 ¹¹ VP	10 Preps	VPK-099
ViraBind™ Adenovirus Purification Kit	2.5 x 10 ¹² VP	10 Preps	VPK-100

QuickTiter™ Adenoviral Titer & Quantitation Kits

Accurate measurement of virus titer is critical for viral gene delivery. Traditional plaque-forming unit (PFU) assays are long and have high inter-assay variability. The QuickTiter™ Adenovirus Titer Kits provide a complete system to functionally titer virus infectivity with greater accuracy in a fraction of the time. The assays may be used with any adenovirus system that can amplify in 293 cells. Assays are available for ICC staining or 96-well ELISA.

For a quick test of physical titer, our QuickTiter™ Adenovirus Quantitation Kit measures the concentration of your adenovirus prep in about one hour.



QuickTiter™ Adenovirus Titer Immunoassay Kit. 293AD cells (p. 42) were infected with different dilutions of purified Ad-β-Gal for 48 hours. Immunostaining was performed according to the assay protocol. X-gal staining was performed with β-Galactosidase Staining Kit (p. 106).

- **Faster, More Accurate and Precise:** Compared to traditional plaque-forming unit assays
- **User-Friendly Protocol:** No agar overlay steps
- **Versatile:** Recognize all 41 adenovirus serotypes

Recent Product Citations

1. Nuche-Berenguer, B. et al. (2015). Elucidation of the roles of the Src kinases in pancreatic acinar cell signaling. *J. Cell Biochem.* **116**:22-36. (VPK-106)
2. Sanchez-Lugo, Y.E. et al. (2014). CXCL10/XCL1 fusokine elicits in vitro and in vivo chemotaxis. *Biotechnol. Lett.* 10.1007/s10529-014-1746-1. (VPK-106)
3. Gibot, L. et al. (2015). Cell-based approach for 3D reconstruction of lymphatic capillaries in vitro reveals distinct functions of HGF and VEGF-C in lymphangiogenesis. *Biomaterials* 10.1016/j.biomaterials.2015.11.027. (VPK-109)
4. Herath, S. et al. (2015). Analysis of T cell responses to chimpanzee adenovirus vectors encoding HIV gag-pol-nef antigen. *Vaccine* 10.1016/j.vaccine.2015.10.111. (VPK-109)
5. Yang, Y. et al. (2015). RGD-modified oncolytic adenovirus exhibited potent cytotoxic effect on CAR-negative bladder cancer-initiating cells. *Cell Death Dis.* **6**:e1760. (VPK-109)
6. Nakao, S. et al. (2015). Stimulus-dependent regulation of nuclear Ca²⁺ signaling in cardiomyocytes: a role of neuronal calcium sensor-1. *PLoS One* **10**:e0125050. (VPK-109)
7. Patsouris, D. et al. (2014). Insulin resistance is associated with MCP1-mediated macrophage accumulation in skeletal muscle in mice and humans. *PLoS One* **9**:e110653. (VPK-109)
8. Barrett, A. et al. (2014). A crucial role for DOK1 in PDGF-BB-stimulated glioma cell invasion through p130Cas and Rap1 signalling. *J. Cell Sci.* **127**:2647-2658. (VPK-109)
9. Garcia-Pascual, C.M. et al. (2015). Evaluation of the potential therapeutic effects of a double-stranded RNA mimic complexed with polycations in an experimental mouse model of endometriosis. *Fertil. Steril.* 10.1016/j.fertnstert.2015.07.1147. (VPK-110)
10. Gibson, H. et al. (2015). Immunotherapeutic intervention with oncolytic adenovirus in mouse mammary tumors. *Oncol Immunology* **4**:e984523. (VPK-110)

Selection Guide for QuickTiter™ Adenoviral Quantitation Kits

	QuickTiter™ Adenovirus Titer Immunoassay Kit	QuickTiter™ Adenovirus Titer ELISA Kit	QuickTiter™ Adenovirus Quantitation Kit
Functional or Physical Titer	Functional (Infectious units)	Functional (Infectious units)	Physical (Viral particles)
Assay Time	2.5 days	2.5 days	45-60 minutes
Assay Principle	Antibody-based	Antibody-based	Total nucleic acid content
Detection Method	Immunocytochemical staining	Colorimetric (ELISA) plate reader	Fluorescence plate reader
Key Benefit	Accuracy	Accuracy	Speed

Product Name	Detection	Size	Catalog Number
QuickTiter™ Adenovirus Titer Immunoassay Kit	ICC Staining	100 Assays	VPK-109
QuickTiter™ Adenovirus Titer ELISA Kit	Colorimetric	2 x 96 Assays	VPK-110
QuickTiter™ Adenovirus Quantitation Kit	Fluorometric	20 Assays	VPK-106

Rapid Replication Competent Adenovirus (RCA) Assay Kit

Traditionally RCA (replication competent adenovirus) is measured in permissive cells by a plaque-forming unit (PFU) assay which takes 10-14 days. This kit uses the assay principles of the QuickTiter™ Adenovirus Titer Immunoassay Kit (see [page 55](#)), but is designed specifically to measure the level of replication-competent virus in your adenoviral prep.

- **Faster Results:** 2.5 days vs. 10 days with the traditional plaque assay
- **Versatile:** Recognizes all 41 adenovirus serotypes

Recent Product Citation

Dubuisson, O. et al. (2015). Accurate identification of neutralizing antibodies to adenovirus Ad36, a putative contributor of obesity in humans. *J. Diabetes Complications* 29:83-87.

Product Name	Detection	Size	Catalog Number
Rapid RCA Assay Kit	ICC Staining	30 Assays	VPK-111
		5 x 30 Assays	VPK-111-5

ViraDuctin™ Adenovirus Transduction Reagent, CAR-Independent



Adenovirus infection of target cells is mediated largely by the coxsackievirus-adenovirus receptor (CAR). Generally adenoviral transduction of many immortalized cell lines proceeds with a high level of efficiency. However, in many primary cells this receptor is either absent or present at extremely low-levels. This can reduce the efficiency of adenovirus transduction into your cell of choice.

ViraDuctin™ Adenovirus Transduction Reagent is designed specifically to increase the efficiency of adenoviral transduction, without regard to the level of CAR expression on the surface of the target cells.

- **Higher Transduction Efficiency:** Up to 12-fold increase in adenoviral uptake
- **User-Friendly:** Short incubation step prior to host cell infection
- **Versatile:** Ideal for target cells expressing little or no CAR, but may also improve transduction efficiency for CAR-expressing cells

Recent Product Citations

1. Haidari, M. et al. (2014). Disruption of endothelial adherens junctions by high glucose is mediated by protein kinase C-β-dependent vascular endothelial cadherin tyrosine phosphorylation. *Cardiovasc. Diabetol.* 13:112.
2. Haidari, M. et al. (2012). Integrin α2β1 mediates tyrosine phosphorylation of vascular endothelial cadherin induced by invasive breast cancer cells. *J. Biol. Chem.* 287:32981-32992.
3. Ackerman, W. et al (2008). Nuclear Factor-kappa B regulates inducible prostaglandin E synthase expression in human amnion mesenchymal cells. *Biol. Reprod.* 78:68-76.



Enhanced Transduction using ViraDuctin™ Adenovirus Transduction Reagent. Infection of NIH3T3 cells with recombinant Ad-β-gal (ADV-002). **Top:** X-gal staining under microscope. **Bottom:** scoring of infection with ViraDuctin™ reagent as a percentage of infection with control.

Product Name	Size*	Catalog Number
ViraDuctin™ Adenovirus Transduction Reagent (CAR-Independent)	10 Transductions	AD-200
	50 Transductions	AD-201

*Based on using 6-well plates or 35mm culture dishes; may also be used with 96-,24- or 12-well plates or 60mm or 100mm dishes.

Lentiviral Expression Kits & Reagents

As a sub-class of retroviruses, lentiviruses based on HIV-1 have the unique advantage of being able to infect both proliferating and non-proliferating cells, and they can be used for both transient and stable gene expression.

We offer a complete workflow solution to your lentiviral expression studies:

- Expression Systems & Vectors
- Premade Controls
- Viral Host Cell Line
- Purification Kits
- Quantitation / Titer Kits
- Transduction Reagents

ViraSafe™ Lentiviral Expression Systems

Lentiviruses based on HIV-1 may infect both dividing and non-dividing cells. Recently developed third-generation lentiviral expression systems have reduced the risk of creating replication-competent virus upon recombination, but the risk is still present.


Our ViraSafe™ Lentiviral Expression Systems provide a safer and more flexible method to package your lentivirus, even compared to other third-generation lentivirus systems.

- **Safer:** 80-90% less sequence homology compared to other 3rd-generation lentiviral systems; ecotropic systems provide even more safety*
- **High Titters:** Incorporates elements that provide titers comparable to other 3rd-generation systems
- **Flexible:** Packaging vectors provided separately for increased safety and optimization of vector ratios

*Lentiviruses made with a ViraSafe™ Ecotropic Expression System will only readily infect mouse and rat cells. Pantropic lentiviruses are VSVG-pseudotyped and may infect cells of any species.

ViraSafe™ Lentiviral Technology is available in three formats (see next two pages for ordering information):

- **Complete Expression Systems:** Include 3 packaging plasmids, expression vector and control vector
- **Packaging Systems:** Include the 3 individual packaging plasmids; ideal if you already have a 3rd-generation lentiviral expression construct
- **Expression Vectors:** 11 cloning vectors to choose from; compatible with any 2nd or 3rd generation packaging system, but produce the highest titers with the ViraSafe™ packaging system



Lentivirus Production using the ViraSafe™ Lentiviral Expression System.

ViraSafe™ Lentiviral Expression Systems, continued

Complete ViraSafe™ Expression Systems include three individual packaging plasmids, an expression vector, and a control vector. Choose an ecotropic system for infection of mouse or rat cells, or a pan-tropic system to produce VSVG-pseudotyped lentivirus for infection of cells from any species.

Recent Product Citation

Davis, M. et al. (2013). RAC1P29S is a spontaneously activating cancer-associated GTPase. *PNAS* **110**:912-917. (VPK-214-PAN)

Product Name	Envelope	Size	Catalog Number
ViraSafe™ Universal Lentiviral Expression System (Promoterless)	Ecotropic	1 Kit	VPK-211-ECO
	Pantropic (VSVG)	1 Kit	VPK-211-PAN
ViraSafe™ Lentiviral Expression System (Puro)	Ecotropic	1 Kit	VPK-212-ECO
	Pantropic (VSVG)	1 Kit	VPK-212-PAN
ViraSafe™ Lentiviral Expression System (Neo)	Ecotropic	1 Kit	VPK-213-ECO
	Pantropic (VSVG)	1 Kit	VPK-213-PAN
ViraSafe™ Lentiviral Expression System (Hygro)	Ecotropic	1 Kit	VPK-214-ECO
	Pantropic (VSVG)	1 Kit	VPK-214-PAN
ViraSafe™ Lentiviral Bicistronic Expression System (Puro)	Ecotropic	1 Kit	VPK-215-ECO
	Pantropic (VSVG)	1 Kit	VPK-215-PAN
ViraSafe™ Lentiviral Bicistronic Expression System (Neo)	Ecotropic	1 Kit	VPK-216-ECO
	Pantropic (VSVG)	1 Kit	VPK-216-PAN
ViraSafe™ Lentiviral Bicistronic Expression System (Hygro)	Ecotropic	1 Kit	VPK-217-ECO
	Pantropic (VSVG)	1 Kit	VPK-217-PAN
ViraSafe™ Lentiviral Bicistronic Expression System (Blasticidin)	Ecotropic	1 Kit	VPK-219-ECO
	Pantropic (VSVG)	1 Kit	VPK-219-PAN
ViraSafe™ shRNA Lentiviral Expression System (Puro)	Ecotropic	1 Kit	VPK-221-ECO
	Pantropic (VSVG)	1 Kit	VPK-221-PAN
ViraSafe™ shRNA Lentiviral Expression System (GFP)	Ecotropic	1 Kit	VPK-222-ECO
	Pantropic (VSVG)	1 Kit	VPK-222-PAN

ViraSafe™ Lentiviral Packaging Systems

ViraSafe™ Packaging Systems contain 3 packaging plasmids for use with any 3rd-generation lentiviral expression vector. These systems are perfect if you already have a lentiviral construct containing your gene of interest.

Recent Product Citation

Vogt, J. et al. (2014). Protein associated with SMAD1 (PAWS1/FAM83G) is a substrate for type I bone morphogenetic protein receptors and modulates bone morphogenetic protein signaling. *Open Bio*. **4**:130210. (VPK-206)

Product Name	Envelope	Size	Catalog Number
ViraSafe™ Lentiviral Packaging System	Ecotropic	1 Kit	VPK-205
	Pantropic (VSVG)	1 Kit	VPK-206

ViraSafe™ Lentiviral Expression Vectors

These lentiviral expression vectors may be used with any 2nd or 3rd generation lentiviral packaging system, but best results are achieved when used with our ViraSafe™ Lentiviral Packaging Systems.

Recent Product Citations

1. Wang, Z. et al. (2015). Elabela-apelin receptor signaling pathway is functional in mammalian systems. *Sci. Rep.* **5**:8170. (VPK-211)
2. Meyer, A.E. (2014). Role of TGF- β receptor III localization in polarity and breast cancer progression. *Mol. Biol. Cell* **25**:2291. (VPK-213)
3. Davis, M. et al. (2013). RAC1P29S is a spontaneously activating cancer-associated GTPase. *PNAS* **110**:912-917. (VPK-214)
4. Belin, B. et al. (2013). Visualization of actin filaments and monomers in somatic cell nuclei. *Mol. Biol. Cell* **24**:982-994. (VPK-219)

Product Name	Cloning Capacity	Size	Catalog Number
pSMPUW Universal Lentiviral Expression Vector (Promoterless)	9.4 kb	10 μ g	VPK-211
pSMPUW-Puro Lentiviral Expression Vector	7.9 kb	10 μ g	VPK-212
pSMPUW-Neo Lentiviral Expression Vector	7.7 kb	10 μ g	VPK-213
pSMPUW-Hygro Lentiviral Expression Vector	7.4 kb	10 μ g	VPK-214
pSMPUW-IRES-Puro Lentiviral Expression Vector	7.8 kb	10 μ g	VPK-215
pSMPUW-IRES-Neo Lentiviral Expression Vector	7.5 kb	10 μ g	VPK-216
pSMPUW-IRES-Hygro Lentiviral Expression Vector	7.3 kb	10 μ g	VPK-217
pSMPUW-IRES-Bsd Lentiviral Expression Vector	7.9 kb	10 μ g	VPK-219
pSMPUW-U6-Puro Lentiviral Expression Vector	7.7 kb	10 μ g	VPK-221
pSMPUW-U6-GFP Lentiviral Expression Vector	7.6 kb	10 μ g	VPK-222

Lentiviral Control Plasmids

Product Name	Size	Catalog Number
pLenti-GFP Lentiviral Control Vector	10 μ g	LTV-400
pSMPUW-GFP-Puro Lentiviral Control Vector	10 μ g	LTV-401
pSMPUW-MNDnLacZ Lentiviral Control Vector	10 μ g	LTV-402
pLenti-RFP-Puro Lentiviral Control Vector	100 μ L	LTV-403

Premade Reporter Lentivirus Controls

Product Name	Concentration	Size	Catalog Number
GFP Lentivirus Control	1 x 10 ⁶ TU/mL	200 μ L	LTV-300
RFP Lentivirus Control	1 x 10 ⁶ TU/mL	200 μ L	LTV-301

293LTV Lentiviral Cell Line

Our 293LTV cell line was selected from the parental 293T cell line for firmer attachment to culture plates and larger, rounder morphology for greater lentiviral production.

Recent Product Citations

1. Billcliff, P.G. et al. (2015). OCRL1 engages with the F-BAR protein pacsin 2 to promote biogenesis of membrane trafficking intermediates. *Mol. Biol. Cell* **10**.1091/mbc.E15-06-0329.
2. Latta, C.H. et al. (2015). Determining the role of IL-4 induced neuroinflammation in microglial activity and amyloid- β using BV2 microglial cells and APP/PS1 transgenic mice. *J. Neuroinflamm.* **12**:41.

Product Name	Size	Catalog Number
293LTV Cell Line	1 x 10 ⁶ Cells	LTV-100


QuickTiter™ Lentivirus Titer / Quantitation Kits

Measuring lentiviral titer is important prior to infection of your target cells, and one of the most published methods is the p24 ELISA. Our traditional p24 ELISA kit provides a quick, convenient way to quantify the concentration of your HIV-1 based lentivirus.

One disadvantage of using a traditional p24 ELISA to quantify lentivirus is the overexpression of p24 during lentiviral packaging. Free p24 protein may account for a substantial portion of total p24 in lentiviral supernatant. The traditional p24 ELISA detects both virus-associated p24 and free p24 generated by 293T cells during transient transfection. Our QuickTiter™ Lentivirus Titer Kit minimizes the overestimation of p24 in lentivirus supernatant. Our proprietary technology separates the lentivirus-associated p24 from free p24 protein prior to performing the ELISA.

If you need a very quick estimate of your lentiviral concentration, try the QuickTiter™ Lentivirus Quantitation Kit. This kit specifically measures the viral nucleic acid content of purified virus or unpurified viral supernatant. This method is ideal for a quick measurement of viral titer, either before or after purification of your lentivirus.

- **More Accurate:** Exclusive technology in QuickTiter™ Lentivirus Titer Kit minimizes overestimation of virus titer
- **User-Friendly:** Read results on a standard microplate reader




Assay Principle for the QuickTiter™ Lentivirus Titer Kit. Lentivirus particles are packaged with p24 protein, but additional free p24 protein is present in viral supernatant. A traditional p24 ELISA detects both sources of p24 which overestimates viral titer. The QuickTiter™ Lentivirus Titer Kit uses technology to pull the virus out of solution prior to quantitation for a more accurate viral titer.


Selection Guide for QuickTiter™ Lentivirus Quantitation & Titer Kits

	QuickTiter™ Lentivirus Titer Kit (Lentivirus-Associated p24 ELISA)	QuickTiter™ Lentivirus Quantitation Kits (Traditional p24 ELISA)	QuickTiter™ Lentivirus Quantitation Kit
Assay Principle	p24 ELISA with proprietary technology to separate free p24 from viral p24	p24 ELISA	Measures nucleic acid content
Suitable Viruses	Recombinant HIV-1	Recombinant or native HIV-1	HIV-1, FIV, SIV
Detection Method	Colorimetric (ELISA) plate reader	Colorimetric (ELISA) plate reader	Fluorescence plate reader
Key Benefit	Accuracy	Most Published	Speed (45-60 min.)

QuickTiter™ Lentivirus Titer / Quantitation Kits, continued



Free p24 Does not Complex with ViraBind™ Reagents. Recombinant p24 was diluted in culture medium and treated with Vira-Bind™ Lentivirus Reagents A and B found in the QuickTiter™ Lentivirus Titer Kit. The amount of p24 in the supernatant and the pellet was measured according to the assay protocol.



p24 Titer of GFP Lentiviral Supernatant. GFP lentiviral construct was cotransfected with a packaging mix into 293 cells. The conditioned medium was harvested 48 hrs after transfection and used to further infect 293 cells. The p24 level of the diluted lentiviral supernatant (1:10 dilution) was determined as described in the assay protocol.

Recent Product Citations

- Loperfido, M. et al. (2015). piggyBac transposons expressing full-length human dystrophin enable genetic correction of dystrophic mesoangioblasts. *Nucleic Acids Res.* 10.1093/nar/gkv1464. (VPK-107)
- Chen, P.Y. et al. (2015). Endothelial-to-mesenchymal transition drives atherosclerosis progression. *J. Clin. Invest.* 10.1172/JCI82719. (VPK-107)
- Feng, Y. et al. (2015). Natural polymorphisms and oligomerization of human APOBEC3H contribute to single-stranded DNA scanning ability. *J. Biol. Chem.* 10.1074/jbc.M115.666065. (VPK-107)
- Noh, K.M. et al. (2015). ATX tolerates activity dependent histone H3 methyl/phos switching to maintain repetitive element silencing in neurons. *PNAS* 112:6820-6827. (VPK-107)
- Chaubaud, M. et al. (2015). Cell migration and antigen capture are antagonistic processes coupled by myosin II in dendritic cells. *Nat. Commun.* 6:7526. (VPK-107)
- Li, L. et al. (2015). Mammalian target of rapamycin overexpression antagonizes chronic hypoxia-triggered pulmonary arterial hypertension via the autophagic pathway. *Int. J. Mol. Med.* 36:316-322. (VPK-107)
- Zhao, S. et al. (2015). The DEAD-box RNA helicase 5 positively regulates the replication of porcine reproductive and respiratory syndrome virus by interacting with viral Nsp9 in vitro. *Virus Res.* 195:217-224. (VPK-107)
- Haggani, A.A. et al. (2015). Central memory DE4+ cells are preferential targets of double infection by HIV-1. *Viral. J.* 12:184. (VPK-108-H)
- Oh, S.M. et al. (2015). Combined Nurr1 and Foxa2 roles in the therapy of Parkinson's disease. *EMBO Mol. Med.* 10.15252/emmm.201404610. (VPK-108-H)
- Tilton, C.A. et al. (2015). A combination HIV reporter virus system for measuring post-entry event efficiency and viral outcome in primary CD4+ T cell subsets. *J. Virol. Methods* 195:164-169. (VPK-108-H)
- Lucera, M.B. et al. (2014). The histone deacetylase inhibitor Vorinostat (SAHA) increases the susceptibility of uninfected CD4+ T cells to HIV by increasing the kinetics and efficiency of postentry viral events. *J. Virol.* 88:10803-10812. (VPK-108-H)
- Lambert, M.P. et al. (2015). Intramedullary megakaryocytes internalize released platelet factor 4 (PF4) and store it in alpha granules. *J. Thromb. Haemost.* 10.1111/jth.13069. (VPK-112)
- Yadavilli, S. et al. (2015). The emerging role of NG2 in pediatric diffuse intrinsic pontine glioma. *Oncotarget* 6:12141-12155. (VPK-112)
- Ahmad, M.A. et al. (2015). Label-free capacitance-based identification of viruses. *Sci. Rep.* 5:9809. (VPK-112)
- Tang, X. et al. (2015). The advantages of PD1 activating chimeric receptor (PD1-ACR) engineered lymphocytes for PDL1+ cancer therapy. *Am. J. Transl. Res.* 7:460-473. (VPK-112)
- Shin, H.S. et al. (2015). Crosstalk among IL-23 and DNAX activating protein of 2 kDa-dependent pathways promotes osteoclastogenesis. *J. Immunol.* 194:316-324. (VPK-112)

Product Name	Detection	Size	Catalog Number
QuickTiter™ Lentivirus Titer Kit (Lentivirus-Associated HIV p24 ELISA)	Colorimetric	96 Assays	VPK-107
		5 x 96 Assays	VPK-107-5
QuickTiter™ Lentivirus Quantitation Kit (HIV-1 p24 ELISA)	Colorimetric	96 Assays	VPK-108-H
		5 x 96 Assays	VPK-108-H-5
QuickTiter™ Lentivirus Quantitation Kit	Fluorometric	20 Assays	VPK-112

ViraBind™ Lentivirus Purification Kit

Ultracentrifugation methods used for lentiviral supernatants are tedious and time-consuming and usually only partially purify your virus.

The ViraBind™ Lentivirus Purification Kit produces purified lentivirus with extremely high titer without the need for ultracentrifugation. The kit uses a proprietary syringe filter for highly pure lentivirus preps.

- **High Viral Yield:** >90% recovery
- **High Quality:** Provides quality of CsCl procedures,
- **Faster Results:** Quick 2 hour protocol

Recent Product Citation

Zhou, C. et al. (2015). Lhx8 mediated Wnt and TGFβ pathways in tooth development and regeneration. *Biomaterials* 10.1016/j.biomaterials.2015.06.004.

Product Name	Capacity/Prep	Size	Catalog Number
ViraBind™ Lentivirus Purification Kit	2 x 10 ⁸ IFU	10 Preps	VPK-104

ViraDuctin™ Lentivirus Transduction Kit


Lentivirus transduction efficiency is typically low. Additives such as Polybrene® can boost transduction efficiencies, but even then only a small fraction of lentiviral vectors can transduce many target cell lines.

Our ViraDuctin™ Lentivirus Transduction Kit provides superior transduction efficiencies in a variety of cell lines, even when compared to transductions in the presence of Polybrene®.

Recent Product Citations

1. Osorio, L.A. et L. (2015). SNAIL transcription factor increases the motility and invasive capacity of prostate cancer cells. *Mol. Med. Rep.* 13:778-786.
2. Kandasamy, K. et al. (2015). Changes in endothelial Cx43 expression inversely correlates with microvessel permeability and VE-cadherin expression in endotoxin challenged lungs. *Am. J. Physiol. Lung Cell Mol. Physiol.* 10.1152/ajplung.00211.2014.
3. Abel, E.V. et al. (2014). The Notch pathway is important in maintaining the cancer stem cell population in pancreatic cancer. *PLoS One* 9:e91983.
4. Ozelo, M.C. et al. (2014). Omental implantation of BOECs in hemophilia dogs results in circulating FVIII antigen and a complex immune response. *Blood* 123:4045-4053.
5. Rossello, R.A. et al. (2013). Mammalian genes induce partially reprogrammed pluripotent stem cells in non-mammalian vertebrate and invertebrate species. *eLife Sci.* 2:e00036.
6. McEachron, T.A. et al. (2010). Protease-activated receptors mediate crosstalk between coagulation and fibrinolysis. *Blood* 116:5037-5044.
7. Zemskova, M. et al. (2010). p53-dependent induction of prostate cancer cell senescence by the PIM1 protein kinase. *Mol. Cancer Res.* 8:1126-1141.

- **Higher Transduction Efficiency:** 2-6x higher in many cell lines compared to Polybrene
- **More Robust:** Useful for transduction of nonpermissive cells, including primary cells and stem cells



Transduction Efficiencies in Various Cell Lines. NIH3T3 cells, HeLa cells, our own 293AD cells (page 36) and HT-1080 cells were each seeded at 50,000 cells/well in a 24-well plate overnight. Cells were infected with GFP lentivirus for 48 hours in the presence of Polybrene® or ViraDuctin™ Lentivirus Transduction Kit. For each cell line, fluorescence levels using the ViraDuctin™ system are depicted relative to a normalized fluorescence of 100 for Polybrene®.

Polybrene is a registered trademark of Abbott Laboratories.

Product Name	Size*	Catalog Number
ViraDuctin™ Lentivirus Transduction Kit	40 Transductions	LTV-200
	200 Transductions	LTV-201

*Based on a 24-well plate. Can also be used with 96-well, 12-well or 6-well plates, as well as 60mm or 100mm dishes. See product insert.

Retroviral Expression Kits & Reagents

Traditional retroviral vectors based on MMLV are useful for integrating genetic material into the host cell genome. However, retrovirus titer tends to be significantly lower than that of adenovirus, which can lead to a lower infection efficiency.

Our retroviral reagents and kits incorporate technologies that increase your chances of successful retroviral expression. We offer a comprehensive solution from start to finish:

- Retroviral Expression Systems
- Retroviral Packaging Cell Lines
- Retroviral Cloning & Expression Vectors
- Gene-Specific Retroviral Vectors
- Concentration / Purification Kits
- Quantitation / Titer Kits
- Transduction Reagents

Platinum Retroviral Packaging Cell Lines

Generate high titers of recombinant retrovirus with a single plasmid transfection* using these extremely powerful, stable cell lines. Platinum Retroviral Packaging Cells are based on the 293T cell line and exhibit greater stability and produce higher yields of retroviral structure proteins, resulting in higher retroviral titers.

The Platinum cell lines were invented in the laboratory of Dr. Toshio Kitamura at the University of Tokyo and are available exclusively from Cell Biolabs. They were first described in the following paper:

Morita, S. et al. (2000). *Gene Therapy* 7:1063-1066.

Recent Product Citations

- Jani, R.A. et al. (2015). STX13 regulates cargo delivery from recycling endosomes during melanosome biogenesis. *J. Cell Sci.* 128:3263-3276. (RV-101)
- Song, A. et al. (2015). Molecular changes associated with acquired resistance to crizotinib in ros1-rearranged non-small cell lung cancer. *Clin. Cancer Res.* 21:2379-2387. (RV-101)
- Ogi, H. et al. (2015). The oncogenic role of the cochaperone Sgt1. *Oncogenesis* 4:e149. (RV-101)
- Chinyenetere, F. et al. (2015). Mice null for the deubiquitinase USP18 spontaneously develop leiomyosarcomas. *BMC Cancer* 15:886. (RV-102)
- Fuerstenau-Sharp, M. et al. (2015). Generation of highly purified human cardiomyocytes from peripheral blood mononuclear cell-derived induced pluripotent stem cells. *PLoS One* 10:e0126596. (RV-102)
- Kuroda, M. et al. (2015). Interaction between TIM-1 and NPC1 is important for cellular entry of Ebola virus. *J. Virol.* 89:6481-6493. (RV-103)

Not sure which Platinum Expression System is right for you? See the table below for a selection guide based on the host species of your target cell.

	Plat-A Cells (Amphotropic)	Plat-E Cells (Ecotropic)	Plat-GP Cells (Pantropic*)
Human	+++	N.S.	+++
Mouse	+++	+++	+++
Rat	+++	+++	+++
Monkey	+++	N.S.	+++
Cat	+++	N.S.	+++
Dog	+++	N.S.	+++
Hamster	+	N.S.	+++
Bird	N.S.	N.S.	+++
Fish	N.S.	N.S.	+++
Frog	N.S.	N.S.	+++
Insect	N.S.	N.S.	+++
Mollusk	N.S.	N.S.	+++

*Plat-GP cells must be co-transfected with a pantropic envelope protein such as VSV-G.
N.S. = Not Suitable

Suitability of Platinum Retroviral Packaging Cell Lines by Host Species.

Product Name	Size	Catalog Number
Platinum-E Retroviral Packaging Cell Line, Ecotropic	$\geq 3 \times 10^6$ cells	RV-101
Platinum-A Retroviral Packaging Cell Line, Amphotropic	$\geq 3 \times 10^6$ cells	RV-102
Platinum-GP Retroviral Packaging Cell Line, Pantropic	$\geq 3 \times 10^6$ cells	RV-103
pVSV-G Packaging Vector	10 μ g	RV-110

Platinum Retroviral Packaging Cells and Expression Systems

Our Platinum Retroviral Expression Systems incorporate superior packaging cell lines and vector technologies to produce high-titer virus with a single plasmid transfection. Each Platinum Expression System includes one of our exclusive Platinum Packaging Cell Lines which stably express the gag and pol genes. In the Ecotropic and Amphotropic systems, the packaging cells also express the envelope protein.* Simply clone your gene of interest into the vector provided and transfect into the Platinum cells.


- **Higher Viral Yields:** Average titer 10^7 infectious units/mL with transient transfection
- **Longer Stability:** Expression up to 4 months in the presence of drug selection
- **Optimized Systems:** 3 packaging cell lines for infection of various species; 3 vector backbones (two specifically for infection of stem cells)
- **Flexible:** Order complete systems or cells and vectors separately

*Pantropic systems require co-transfection with the provided VSVG envelope vector.

Recent Product Citations

1. Aoi, N. et al. (2012). 1 α ,25-dihydroxyvitamin D3 modulates the hair-inductive capacity of dermal papilla cells: therapeutic potential for hair regeneration. *Stem Cells Trans Med.* 1:615-626. (VPK-301)
2. Saito, Y. et al. (2015). Enhancement of spontaneous activity by HCN4 overexpression in mouse embryonic stem cell-derived cardiomyocytes—a possible biological pacemaker. *PLoS One* 10:e0138193. (VPK-302)
3. Yamaguchi, J. et al. (2015). Inflammation and hypoxia linked to renal injury by CCATT/enhancer-binding protein δ . *Kidney Int.* 10.1038/ki.2015.21. (VPK-302)
4. Wang, N. et al. (2013). Lacritin rescues stressed epithelia via rapid forkhead box O3 (FOXO3)-associated autophagy that restores metabolism. *J. Biol. Chem.* 288:18146-18161. (VPK-302)
5. Kishida, T. et al. (2015). Reprogrammed functional brown adipocytes ameliorate insulin resistance and dyslipidemia in diet-induced obesity and type 2 diabetes. *Stem Cell Reports* 10.1016/j.stemcr.2015.08.007. (VPK-303, VPK-305)

Each system contains a packaging cell line, an expression vector, and GFP control vector. Our pan-tropic systems also contain a VSVG envelope vector.



Retrovirus Production Using the Platinum Expression Systems (Ecotropic and Amphotropic).

Product Name	Expression Vector	Packaging Cell	Catalog Number
Platinum Retroviral Expression System, Ecotropic	pMXs-Puro	Plat-E	VPK-300
Platinum Retroviral Expression System, Amphotropic	pMXs-Puro	Plat-A	VPK-301
Platinum Retroviral Expression System, Pantropic	pMXs-Puro	Plat-GP	VPK-302
Platinum ES/EC Retroviral Expression System, Ecotropic	pMCs-Puro	Plat-E	VPK-303
Platinum ES/EC Retroviral Expression System, Amphotropic	pMCs-Puro	Plat-A	VPK-304
Platinum ES/EC Retroviral Expression System, Pantropic	pMCs-Puro	Plat-GP	VPK-305
Platinum HSC Retroviral Expression System, Ecotropic	pMYs-Puro	Plat-E	VPK-306
Platinum HSC Retroviral Expression System, Amphotropic	pMYs-Puro	Plat-A	VPK-307
Platinum HSC Retroviral Expression System, Pantropic	pMYs-Puro	Plat-GP	VPK-308

Retroviral Cloning & Expression Vectors

Our Retroviral Expression Vectors are based on backbones derived from Moloney murine leukemia virus (MMLV). We offer the traditional pBABE system and the novel pMXs system, which has been shown to be useful in induced pluripotent stem cell (iPS) studies. pMYs vectors are optimal for use with hematopoietic stem cells, and pMCs vectors are optimal for ES and EC cells. All cloning vectors are supplied as 10 µg in TE buffer.

Retroviral Cloning Vectors for General Gene Expression (driven by 5' LTR)

Recent Product Citations

1. Diep, C.M. et al. (2015). Retroviral expression of human cystatin genes in HeLa cells. *Methods Mol. Bio.* **1249**:121-131. (RTV-001-PURO)
2. Nakamura, H. et al. (2015). Genomic spectra of biliary tract cancer. *Nat. Genet.* **47**:1003-1010. (RTV-010)
3. Yeon, J.T. et al. (2015). Arginase 1 is a negative regulator of osteoclast differentiation. *Amino Acids* **10.1007/s00726-015-2112-0**. (RTV-012)
4. Kageyama-Yahara, N. et al. (2014). Gli regulates MUC5AC transcription in human gastrointestinal cells. *PLoS One* **9**:e106106. (RTV-013)
5. Koso, H. et al. (2014). Identification of FoxR2 as an oncogene in medulloblastoma. *Cancer Res.* **74**:2351-2361. (RTV-014, RTV-015)

Vector Name	Cloning Capacity	Catalog Number
pBABEhygro	5.6 kb	RTV-001-HYGRO
pBABEneo	5.9 kb	RTV-003
pBABEpuro	6 kb	RTV-001-PURO
pBABEzeo	6.3 kb	RTV-004
pMXs	5.4 kb	RTV-010
pMXs-IRES-Bsd	5.6 kb	RTV-016
pMXs-IRES-GFP	5.3 kb	RTV-013
pMXs-IRES-Neo	5.2 kb	RTV-015
pMXs-IRES-Puro	5.4 kb	RTV-014
pMXs-Neo	3.8 kb	RTV-011
pMXs-Puro	4.4 kb	RTV-012
pMZs	5.3 kb	RTV-030

Retroviral Cloning Vector for miRNA

Recent Product Citation

Mansour, M. et al. (2013). The TAL1 complex targets the FBXW7 tumor suppressor by activating miR-223 in human T cell acute lymphoblastic leukemia. *J. Exp. Med.* **210**:1545-1557. (RTV-017)

Vector Name	Cloning Capacity	Catalog Number
pMXs-miR-GFP/Puro	4.2 kb	RTV-017

Retroviral Cloning Vectors with Strong Promoters for Overexpression

Vector Name	Cloning Capacity	Catalog Number
pMXs-CAG	5.2 kb	RTV-064
pMXs-CMV	5.5 kb	RTV-065
pMXs-EF1 α	5.5 kb	RTV-063
pMXs-EF1-Bsd	4.2 kb	RTV-062
pMXs-EF1-GFP	3.9 kb	RTV-061
pMXs-EF1-Puro	4 kb	RTV-060
pMXs-SR α	5.4 kb	RTV-066

Retroviral Cloning Vectors for use with ES/EC Cells

Recent Product Citation

Mochizunki, Y. et al. (2013). Phosphatidylinositol 3-phosphate myotubularin-related protein 6 (MTMR6) is regulated by small GTPase Rab1b in the early secretory and autophagic pathways. *J. Biol. Chem.* **288**:1009-1021. (RTV-041)

Vector Name	Cloning Capacity	Catalog Number
pMCs-IRES-GFP	5.2 kb	RTV-040
pMCs-Puro	4.3 kb	RTV-041

Retroviral Cloning Vectors for use with Hematopoietic Cells

Recent Product Citation

Xiao, X. et al. (2015). GITR subverts Foxp3+ Tregs to boost Th9 immunity through regulation of histone acetylation. *Nat. Commun.* **6**:8266. (RTV-021)

Vector Name	Cloning Capacity	Catalog Number
pMYs	5.2 kb	RTV-020
pMYs-IRES-GFP	5.2 kb	RTV-021
pMYs-IRES-Neo	5.2 kb	RTV-023
pMYs-IRES-Puro	5.4 kb	RTV-022
pMYs-Puro	4.3 kb	RTV-024

Retroviral Cloning Vectors for shRNA

Vector Name	Cloning Capacity	Catalog Number
pMXs-U6-GFP	5 kb	RTV-071
pMXs-U6-Puro	5.1 kb	RTV-070
pMXs-U6-Puro-shGFP		RTV-055
pMXs-U6-Puro-shLuc		RTV-056

Retroviral Packaging Vectors and Cells

Recent Product Citations

- Amagai, Y. et al. (2015). A point mutation in the extracellular domain of KIT promotes tumorigenesis of mast cells via ligand-independent auto-dimerization. *Sci. Rep.* **5**:9775. (RV-110)
- Zhang, T. et al. (2015). Homoharringtonine binds to and increases myosin-9 in myeloid leukemia. *Br. J. Pharmacol.* 10.1111/bph.13359. (RV-110)
- Okamoto, K. et al. (2012). Dengue virus strain DEN2 16681 utilizes a specific glycochain of syndecan-2 proteoglycan as a receptor. *J. Gen. Virol.* **93**:761-770. (RV-110, RV-111)
- Ng, A.J. et al. (2015). The DNA helicase Recq14 is required for normal osteoblast expansion and osteosarcoma formation. *PLoS Genet.* 10:e1005160. (RV-112)

Product Name	Size	Catalog Number
pCMV-10A1 Envelope Vector	100 µL	RV-114
pCMV-Ampho Envelope Vector	100 µL	RV-113
pCMV-Eco Envelope Vector	100 µL	RV-112
pCMV-Gag-Pol Retroviral Vector	10 µg	RV-111
pCMV-VSV-G Envelope Vector	10 µg	RV-110

293RTV Cell Line

Our 293RTV cells are derived from the 293 parental cell line, but are selected for firmer attachment to culture plates, faster growth and higher yields of retrovirus produced.

Product Name	Size	Catalog Number
293RTV Cell Line	≥1 x 10 ⁶ cells	RV-100

Gene-Specific Recombinant Retroviral Vectors

These constructs are based on backbones derived from MMLV. Vectors with GFP or stem cell factors are supplied as 10 µg of plasmid in TE buffer. All other vectors are supplied as 100 µL of bacterial glycerol stock. Product listing continues on the following pages.

Cell Cycle

Recent Product Citation

Huang, J. et al. (2009). Regulation of the leucocyte chemoattractant receptor FPR in glioblastoma cells by cell differentiation. *Carcinogenesis* **30**(2):348-355. (RTV-401)

Target Name	Vector Backbone	Catalog Number
c-Abl	pBABEpuro	RTV-402
c-Abl-TM	pBABEpuro	RTV-403
c-Abl (1-565)	pBABEpuro	RTV-404
c-Abl (1-958)	pBABEpuro	RTV-405
hTERT	pBABEhygro	RTV-007
	pBABEneo	RTV-005
	pBABEpuro	RTV-006
p53	pBABEpuro	RTV-401

Autophagy

This vector is supplied with a separate pMXs-GFP control vector at no additional cost.

Target Name	Vector Backbone	Catalog Number
GFP-LC3	pMXs	RTV-801

Reporter Genes

Recent Product Citations

- Hrdlickova, R. et al. (2012). Alternatively spliced telomerase reverse transcriptase variants lacking telomerase activity stimulate cell proliferation. *Mol. Cell Biol.* **32**:4283-4296. (RTV-002)
- Malaver-Ortega, L.F. et al. (2013). Inducing pluripotency in cattle. *Methods Mol. Biol.* 10.1007/978-1-4939-2848-4_6. (RTV-050)

Target Name	Vector Backbone	Catalog Number
GFP	pBABE	RTV-002
GFP	pMCs	RTV-051
GFP	pMX	RTV-050
GFP	pMYs	RTV-052
GFP-Puro	pMX	RTV-053

Gene-Specific Recombinant Retroviral Vectors, continued

Cytoskeleton Regulation

Recent Product Citation

Zhao, B. et al. (2012). TNF-induced osteoclastogenesis and inflammatory bone resorption are inhibited by transcription factor RBP-J. *J. Exp. Med.* **209**:2467-2483. (RTV-101)

Target Name	Vector Backbone	Mutation State	Catalog Number
Cdc42	pBABEhygro	L61	RTV-203
K-Ras	pBABEpuro	N/A	RTV-220
	pWZLhygro	Q61	RTV-221
myr-Rac1	pBABEpuro	N/A	RTV-201
		V12	RTV-206
Rac1	pBABEhygro	V12	RTV-202
N-Ras	pBABEpuro	K61	RTV-222
Rac3	pBABEhygro	V12	RTV-205
Ras	pBABEpuro	V12	RTV-101
	pBABEpuro	V12C40	RTV-104
	pBABEpuro	V12G37	RTV-103
	pBABEpuro	V12S35	RTV-102
RhoA	pBABEhygro	L63	RTV-204

iPS / Stem Cell Factors

Human iPS Genes

Target Name	Vector Backbone	Catalog Number
4-Vector Set*	pMXs	RTV-701-C
6-Vector Set**	pMXs	RTV-709-C
c-Myc	pMXs	RTV-703
Klf4	pMXs	RTV-704
Lin-28	pMXs	RTV-710
NANOG	pMXs	RTV-709
Oct-3/4	pMXs	RTV-701
Sox2	pMXs	RTV-702
p53 shRNA	pRetro	RTV-410

Mouse iPS Genes

Target Name	Vector Backbone	Catalog Number
4-Vector Set*	pMXs	RTV-705-C
6-Vector Set**	pMXs	RTV-711-C
c-Myc	pMXs	RTV-707
Klf4	pMXs	RTV-708
Lin-28	pMXs	RTV-712
NANOG	pMXs	RTV-711
Oct-3/4	pMXs	RTV-705
Sox2	pMXs	RTV-706
p53 shRNA	pRetro	RTV-400

*4-Vector sets contain individual constructs with the following genes: c-Myc, Klf4, Oct-3/4 and Sox2.

**6-Vector sets contain individual constructs with the following genes: c-Myc, Klf4, Oct-3/4, Sox2, Lin-28 and NANOG.

Proteases and Related Molecules

Target Name	Vector Backbone	Catalog Number
uPA	pBABEpuro	RTV-501
uPAR	pBABEhygro	RTV-502

Recent Product Citation

Gutova, M. et al (2008). Urokinase plasminogen activator and urokinase plasminogen activator receptor mediate human stem cell tropism to malignant solid tumors. *Stem Cells* **26**:1406-1413. (RTV-501, RTV-502)

Gene-Specific Recombinant Retroviral Vectors, continued

MAP Kinase Signaling

Vector Name	Vector Backbone	Mutation State	Catalog Number
ERK2	pBABEhygro	Dominant Negative	RTV-109
JNK1	pBABEpuro	Dominant Negative	RTV-110
MAPKAPK2	pBABEpuro	Constitutively Active	RTV-118
	pBABEpuro	Dominant Negative	RTV-119
MAPKAPK3	pBABEpuro	Constitutively Active	RTV-120
	pBABEpuro	Dominant Negative	RTV-121
MEK1	pBABEhygro	Constitutively Active	RTV-112
	pBABEhygro	Dominant Negative	RTV-111
MKK3	pBABEpuro	Constitutively Active	RTV-114
	pBABEhygro	Dominant Negative	RTV-115
MKK6	pBABEpuro	Constitutively Active	RTV-116
	pBABEhygro	Dominant Negative	RTV-117
myr-Akt1	pWZLneo	Constitutively Active	RTV-125
p38 α	pBABEhygro	Dominant Negative	RTV-105
p38 β	pBABEhygro	Dominant Negative	RTV-106
p38 γ	pBABEhygro	Dominant Negative	RTV-107
p38 δ	pBABEhygro	Wild Type	RTV-128
PI3K p110 α -CAAX	pWZLneo	Constitutively Active	RTV-124
PRAK	pBABEpuro	Constitutively Active	RTV-122
	pBABEpuro	Dominant Negative	RTV-123
Raf1-CAAX	pWZLneo	Constitutively Active	RTV-113

Transcription Regulation

Target Name	Vector Backbone	Catalog Number
AUF1	pBABEpuro	RTV-305
hnRNPA0	pBABEpuro	RTV-310
hnRNP-A2	pBABEpuro	RTV-340
HuB	pBABEpuro	RTV-302
HuC	pBABEpuro	RTV-303
HuD	pBABEpuro	RTV-301
HuR	pBABEpuro	RTV-304
PABP	pBABEpuro	RTV-307
Stat5A	pMXs	RTV-330
Stat5A(1*6)	pMXs	RTV-331

Recent Product Citation



Yu, Y. et al. (2012). Bcl11a is essential for lymphoid development and negatively regulates p53. *J. Exp. Med.* **209**:2467-2483. (RTV-331)

Target Name	Vector Backbone	Catalog Number
Stat5A-IRES-GFP	pMXs	RTV-332
Stat5A(1*6)-IRES-GFP	pMXs	RTV-333
Stat5B	pMXs	RTV-334
Stat5B(1*6)	pMXs	RTV-335
TIA-1	pBABEpuro	RTV-309
TIAR	pBABEpuro	RTV-308
TTP	pBABEpuro	RTV-306

QuickTiter™ Retrovirus Rapid Quantitation Kit

This kit specifically measures the viral nucleic acid content of purified virus or unpurified viral supernatant. This method is ideal for a quick measurement of viral titer, either before or after purification of your retrovirus.


- **Ultra-fast Results:** 45-60 minute procedure
- **Sensitive:** Limit of detection = 1.5×10^9 VP/mL from 2 mL of retroviral supernatant



Retrovirus RNA Standard Curve. The QuickTiter™ Retrovirus RNA Standard was diluted according to the assay protocol. Fluorescence was measured on a SpectraMax Gemini XS Fluorometer.

Recent Product Citations

1. Manian, K.V. et al. (2015). Understanding the molecular basis of heterogeneity in induced pluripotent stem cells. *Cell Reprogram.* 17:427-440.
2. Ito, T. et al. (2012). Stem cell factor programs the mast cell activation phenotype. *J. Immunol.* 188:5428-5437.



Assay Procedure for the QuickTiter™ Retrovirus Quantitation Kit.

Product Name	Detection	Size	Catalog Number
QuickTiter™ Retrovirus Quantitation Kit	Fluorometric	20 Assays	VPK-120

ViraDuctin™ Retrovirus Transduction Kit

The efficiency of retrovirus transduction can be low compared to other viruses. The rate at which retroviral vectors bind to cells is controlled mostly by diffusion. Additionally, the presence of transduction inhibitors such as proteoglycans and glycosaminoglycans in retroviral supernatants can lead to poor gene transfer. Additives such as Polybrene® can boost transduction efficiencies, but they do not eliminate these transduction inhibitors.

Our ViraDuctin™ Retrovirus Transduction Kit provides superior transduction efficiencies even when compared to transductions in the presence of Polybrene®. A proprietary reagent cocktail forms a supercomplex with the retrovirus which is pelleted away from the supernatant, removing detrimental transduction inhibitors that decrease infection efficiency.

- **More Robust:** Removes harmful transduction inhibitors from retroviral supernatant
- **Higher Transduction Efficiencies:** Compared to infections in the presence of Polybrene or no additives
- **Versatile:** Particularly useful for nonpermissive cells including primary cells and stem cells, but may boost transduction rates in a wide variety of cells

Recent Product Citations

1. Gandhi, M. et al. (2012). Homologous chromosomes make contact at the sites of double-strand breaks in genes in somatic G0/G1-phase human cells. *PNAS* **109**:9454-9459.
2. Miyoshi, N. et al. (2010). Defined factors induce reprogramming of gastrointestinal cancer cells. *PNAS* **107**:40-45.

Product Name	Size*	Catalog Number
ViraDuctin™ Retrovirus Transduction Kit	40 Transductions	RV-200
	200 Transductions	RV-201

*Number of transductions shown is based on use in a 24-well plate. This product may also be used with 96-well, 12-well or 6-well plates, as well as 60mm or 100mm dishes. See product insert for specific details.

Polybrene is a registered trademark of Abbott Laboratories.



Biogenuix Medsystems Pvt. Ltd.

412-B, Jyoti Shikhar Building, District Center,
Janakpuri, New Delhi-110058, INDIA

Phone : +91-11-4875-4875

E-mail : contact@biogenuix.com

Web : www.biogenuix.com